



## Not All Decellularization Processing is the Same...

Remove Donor Cells

Anionic  
Detergent

DNA Removal

Nuclease  
Treatment

dCELL<sup>®</sup> technology

Protease  
Inhibitors

Preserve Tissue  
Structures

Terminally  
Sterilized

Sterility Assurance  
Level  $10^{-6}$


DermaPure<sup>®</sup> is created using innovative dCELL<sup>®</sup> Technology, allowing the dermal tissue to maintain its structure and biomechanical properties. The result is a versatile, intact extracellular matrix.

- Reduces cost and waste.
- Allows for fluid egress.
- Covers a large surface area.
- Eliminates time consuming meshing in the OR.

EXPANDS<sup>1</sup> UP TO **180%**



Non-Expanded DermaPure<sup>®</sup> Meshed



Expanded DermaPure<sup>®</sup> Meshed

## DermaPure<sup>®</sup> and DermaPure<sup>®</sup> Meshed provide the following:

- A cover or barrier to the area needing replacement of integumental tissue.
- Soft tissue with supplemental support, reinforcement, and/or protection from mechanical stress.



## The dCELL<sup>®</sup> difference is clear.

| SKU          | SIZE         | PRODUCT DESCRIPTION                                    |
|--------------|--------------|--------------------------------------------------------|
| 010200HD     | 1 cm X 2 cm  | DermaPure Decellularized Dermal Allograft              |
| 020200HD     | 2 cm x 2 cm  | DermaPure Decellularized Dermal Allograft              |
| 020300HD     | 2 cm x 3 cm  | DermaPure Decellularized Dermal Allograft              |
| 030400HD     | 3 cm x 4 cm  | DermaPure Decellularized Dermal Allograft              |
| 040600HD     | 4 cm x 6 cm  | DermaPure Decellularized Dermal Allograft              |
| 040600HDMESH | 4 cm x 6 cm  | DermaPure Meshed Decellularized Dermal Allograft       |
| 071000HD     | 7 cm x 10 cm | DermaPure Decellularized Dermal Allograft              |
| 0071000HDN   | 7 cm x 10 cm | DermaPure Non-Oriented Decellularized Dermal Allograft |
| 07100HDMESH  | 7 cm x 10 cm | DermaPure Meshed Decellularized Dermal Allograft       |
| 091200HD     | 9 cm x 12 cm | DermaPure Decellularized Dermal Allograft              |

TRX BioSurgery  
Customer Service:  
**(833) 567-0767**

For additional information:  
**tissueregenixus.com**

1. Applicable to DermaPure 7cmx10cm Meshed 3:1; Henderson J, Arya R, Gillespie P. Skingraft meshing, over-meshing and cross-meshing. Int J Surg 10 (2012) 547-550.

Tissue Regenix Wound Care, Inc. is doing business as TRX BioSurgery. TRX BioSurgery, Tissue Regenix Wound Care, Inc. and Tissue Regenix Wound Care Limited are part of the Tissue Regenix Group. DermaPure<sup>®</sup> is a registered trademark of Tissue Regenix Wound Care Limited. dCELL<sup>®</sup> Technology is a registered trademark of Tissue Regenix Limited. TR091 rev03