

The Versatile, Decellularized Dermal Allograft


DERMA [PURE]
dCELL[®] technology

Utilized for pelvic organ prolapse and soft tissue anastomosis procedures.

- Decellularized dermal allograft with an intact extracellular matrix.
- Targeted support for pelvic and soft tissue integrity.
- Zero reported adverse events.



**One Allograft. Multiple Surgical Uses.
DermaPure adapts to your reconstructive needs.**



Pelvic Organ Prolapse (POP) Impacts:

- 1 in 2 women who have given birth
- 2 out of 5 women over the age of 60


Behind every statistic is a woman
whose life is impacted by POP.
She deserves innovation she can trust.

Not All Decellularization Processing is the Same...

**dCELL Technology
is the difference.**

Designed for Structural Integrity, Ease of Use, & Surgical Versatility

Ease of Use



Ambient
temperature
storage



Rinse to
remove
protective
glycerol



Tissue tab
to aid in
orientation



Basement
membrane
& dermal
layer



Dermal
layer on
both sides

Available oriented
and non-oriented

For additional information:
tissueregenixus.com

TRX BioSurgery Customer Service:
(833) 567-0767

Intended Use

When implanted provides soft tissue with supplemental support, reinforcement, and / or protection against mechanical stress.

The intended use and the option of oriented or non-oriented make DermaPure ideal for complex reconstructions, including pelvic floor repair and soft tissue anastomosis.

DermaPure Decellularized Dermal Allograft

SKU	SIZE
040600HD	4 cm x 6 cm
040600HDMESH	4 cm x 6 cm (mesh)
071000HD	7 cm x 10 cm
07100HDMESH	7 cm x 10 cm (mesh)
0071000HDN	7 cm x 10 cm (non-oriented)
091200HD	9 cm x 12 cm
162000HD	16 cm x 20 cm